



Iris
Biotech

CLICK CHEMISTRY

Three Generations

→ **Le Click C'est Chic**

Discover building blocks suitable for 1st, 2nd and 3rd generation Click chemistry – for *in vitro* and *in vivo* applications. Choose one and click it as simply as doing a jigsaw puzzle.

pages 1, 2

1st generation Click chemistry:
CuAAC

[page 1](#)

2nd generation Click chemistry:
SPAAC

[page 1 – 2](#)

3rd generation Click chemistry:
IEDDA

[page 3](#)




Version: IF4_3

Click Chemistry

Three Generations

1st Generation Click Chemistry: CuAAC

Alkynes and azides can undergo a Cu(I)-catalyzed azide-alkyne 1,3-dipolar cycloaddition (CuAAC) to afford 1,4-disubstituted 1,2,3-triazoles. Developed by K. Barry Sharpless and Morton Meldal, this type of chemical transformation was quickly dubbed "Click Chemistry". This reaction has found broad applicability in various fields and is as such currently the most widely used conjugation technique.




The copper-catalyzed azide-alkyne cycloaddition affords the 1,4 disubstituted isomer.

Worth noting is the fact that ruthenium is also able to catalyze a 1,3-dipolar cycloaddition between an azide and an alkyne affording the 1,5-disubstituted regioisomer instead.

However, the presence of copper limits *in vivo* applications due to high cell toxicity, undesired oxidation of proteins and the inhibition of luminescence properties of nanocrystals.

2nd Generation Click Chemistry: SPAAC


The so-called strain-promoted azide-alkyne cycloaddition (SPAAC) is widely utilized as copper-free click reaction in a broad range of research fields. However, when using cyclooctynes, their high reactivity comes with the drawback of causing thermal decomposition and/or oligomerization during storage and experimental handling. Furthermore, competitively occurring addition reaction with biological thiols, e.g. cysteine residues in proteins or peptides, diminish their biorthogonality. In contrast, the relatively new developed 4,8-diazacyclononynes (DACNs) possess high thermal and chemical stability while showing high reactivity and increased hydrophilicity as well as high selectivity towards ynophiles.



Schematic illustration of SPAAC.

CliCr®

A superior class of reagents for metal-free click strain promoted conjugation with azides is based on CliCr®, which is based on the small-molecule TMTH-sulfoximine (TMTHSI). The 7-membered CliCr® ring can be conveniently functionalized with a variety of linkers, e.g. via acylation, sulfonylation, N-alkylation, or carbamoylation. Hence, CliCr® reagents can be used in diverse application, for example for the construction of Antibody-Drug Conjugates (ADCs), for *ex vivo* cell modification (e.g. CAR-T), for small molecule-drug conjugates, for oligonucleotide conjugates as well as for diagnostic labelling of a variety of agents.



Chemical structure of the CliCr® base compound and its derivatization possibilities.

Key benefits are its fast reactivity and increased yield, its hydrophilicity compared to other metal-free click reagents combined with its small size as well as its broad compatibility.


CliCr® is provided under an intellectual property license from Cristal Therapeutics. The trademark CliCr® is the property of Cristal Therapeutics. For information on purchasing a license of CliCr® reagents, contact Cristal Therapeutics via Oxfordlaan 55, 6229 EV Maastricht (The Netherlands) or via info@cristaltherapeutics.com.

Product details

RL-4180 CliCr® base compound

TMTH-Sulfoximine


CAS-No.	2408481-82-1
Formula	C ₁₀ H ₁₇ NOS
Mol. weight	199,31 g/mol



RL-4190 CliCr®-beta-Ala-NH₂*TFA

TMTH-sulfoximine beta-alanine amide TFA salt

Formula	C ₁₃ H ₂₂ N ₂ O ₂ S*CF ₃ COOH
Mol. weight	270,39*114,02 g/mol




RL-4200 CliCr®-Suc

TMTH-sulfoximine succinic acid

CAS-No. 2479971-29-2

Formula C₁₄H₂₁NO₄S

Mol. weight 299,39 g/mol



Product details




RL-4330 CliCr®-Osu

TMTH-sulfoximine succinimidyl ester

CAS-No. 2408481-89-8

Formula C₁₅H₂₀N₂O₅S

Mol. weight 340,39 g/mol



References:

- TMTHSI, a superior 7-membered ring alkyne containing reagent for strain-promoted azide-alkyne cycloaddition reactions; J. Weterings, C. J. F. Rijcken, H. Veldhuis, T. Meulemans, D. Hadavi, M. Timmers, M. Honing, H. Ippel, R. M. J. Liskamp; *Chem. Sci.* 2020; **11**: 9011-9016. <https://doi.org/10.1039/d0sc03477k>
- Exploring the Chemical Properties and Medicinal Applications of Tetramethylthiocycloheptyne Sulfoximine Used in Strain-Promoted Azide-Alkyne Cycloaddition Reactions; M. Timmers, A. Kipper, R. Frey, S. Notermans, M. Voievudskyi, C. Wilson, N. Hentzen, M. Ringle, C. Bovino, B. Stump, C. J. F. Rijcken, T. Vermonden, I. Dijkgraaf, R. Liskamp; *Pharmaceuticals* 2023; **16**: 1155. <https://doi.org/10.3390/ph16081155>
- Specific N-terminal attachment of TMTHSI linkers to native peptides and proteins for strain-promoted azide alkyne cycloaddition; M. Timmers, W. Peeters, N. J. Hauwert, C. J. F. Rijcken, T. Vermonden, I. Dijkgraaf, R. M. J. Liskamp; *Chem. Commun.* 2023; <https://doi.org/10.1039/d3cc03397j>




You need more information about the CliCr® technology?

Watch the recording of our online workshop!



3rd Generation Click Chemistry: IEDDA

Tetrazine ligation represents the option for a copper-free, rapid, and fully biorthogonal type of Click chemistry. Mechanistically, this reaction proceeds via an inverse electron-demand Diels-Alder (IEDDA) cycloaddition reaction between a tetrazine and a strained alkene, such as a *trans*-cyclooctene (TCO), cyclopropane or norbornene, followed by a retro-Diels-Alder reaction under elimination of N₂, the latter rendering the reaction irreversible. As additional benefit, this reaction excels at very low concentrations and can be performed in aqueous media.




Common reaction partners for tetrazines.

There are two main types of tetrazines that are widely applied: 6-Methyl-substituted tetrazines (MeTz) and 6-hydrogen-substituted tetrazines (HTz). MeTz exhibit a high stability while still offering faster reaction kinetics with TCO derivatives than any other biorthogonal reaction pairs (approx. 1000 M⁻¹s⁻¹). HTzs show lower stability and less tolerance to harsh reaction conditions but offer even faster reaction kinetics (up to 30000 M⁻¹s⁻¹) for applications like *in vivo* imaging.



For more detailed information on traditional and copper-free Click conjugation, please see our brochure Click Chemistry.


Free Download is available on our website!



Product details

DBCO derivatives**PEG6740 DBCO-PEG(4)-TFP**


Dibenzocyclooctyne-tetra(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenol ester

Formula $C_{37}H_{38}F_4N_2O_8$
Mol. weight 714,7 g/mol


Product details

**RL-2490 DBCO-mal**


Dibenzocyclooctyne-maleimide

CAS-No. 1395786-30-7
Formula $C_{25}H_{21}N_3O_4$
Mol. weight 427,45 g/mol**PEG6785 DBCO-PEG(36)-MAL**


Dibenzocyclooctyne-36(ethylene glycol)-maleimide

Formula $C_{100}H_{170}N_4O_{41}$
Mol. weight 2084,45 g/mol**PEG6765 DBCO-PEG(36)-TFP**


Dibenzocyclooctyne-36(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenol ester

Formula $C_{101}H_{166}F_4N_2O_{40}$
Mol. weight 2124,41 g/mol**PEG7460 DBCO-PEG(24)-OMe**


alpha-Methoxy-24(ethylene glycol)-amido-dibenzocyclooctyne

Formula $C_{68}H_{114}N_2O_{26}$
Mol. weight 1375,63 g/mol**PEG6750 DBCO-PEG(12)-TFP**


Dibenzocyclooctyne-dodeca(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenol ester

Formula $C_{53}H_{70}F_4N_2O_{16}$
Mol. weight 1067,12 g/mol**RL-2421 DBCO-Sulfo-PEG(4)-NH₂**


Dibenzocyclooctyne-tetra(ethylene glycol)amine

CAS-No. 2055198-05-3
Formula $C_{32}H_{42}N_4O_{10}S$
Mol. weight 674,76 g/mol**RL-2420 DBCO-PEG(4)-NH₂*TFA**


Dibenzocyclooctyne-tetra(ethylene glycol)-amine trifluoro acetic acid salt

CAS-No. 1255942-08-5 net
Formula $C_{29}H_{37}N_3O_6^+ * C_2F_3HO_2^-$
Mol. weight 523,62*114,02 g/mol**RL-2550 DBCO-mPEG (20kDa)**


alpha-Dibenzocyclooctyne-omega-methoxy-poly(ethylene glycol)













CAS-No. 2262541-53-5
Mol. weight 20000 Da**RL-2540 DBCO-mPEG (10kDa)**













alpha-Dibenzocyclooctyne-omega-methoxy-poly(ethylene glycol)

CAS-No. 2262541-53-5
Mol. weight 10000 Da**RL-2530 DBCO-mPEG (5kDa)**

alpha-Dibenzocyclooctyne-omega-methoxy-poly(ethylene glycol)

CAS-No. 2262541-53-5
Mol. weight 5000 Da

		Product details	Product details
RL-2500	DBCO-PEG(4)-mal	<p>Dibenzoazacyclooctyne-tetra(ethylene glycol)-maleimide</p> <p>CAS-No. 1480516-75-3</p> <p>Formula $C_{36}H_{42}N_4O_9$</p> <p>Mol. weight 674,74 g/mol</p>	 
RL-2560	DBCO-mPEG (30kDa)	<p>alpha-Dibenzoazacyclooctyne-omega-methoxy-poly(ethylene glycol)</p> <p>CAS-No. 2262541-53-5</p> <p>Mol. weight 30000 Da</p>	 
RL-2440	DBCO-NHS	<p>Dibenzoazacyclooctyne-carboxylic acid succinimidyl ester</p> <p>CAS-No. 1384870-47-6</p> <p>Formula $C_{25}H_{22}N_2O_5$</p> <p>Mol. weight 430,45 g/mol</p>	 
RL-2430	DBCO-COOH	<p>Dibenzoazacyclooctyne-carboxylic acid</p> <p>CAS-No. 1425485-72-8</p> <p>Formula $C_{21}H_{19}NO_3$</p> <p>Mol. weight 333,38 g/mol</p>	 
RL-2510	DBCO-PEG(4)-OH	<p>Dibenzoazacyclooctyne-tetra(ethylene glycol)</p> <p>CAS-No. 1416711-60-8</p> <p>Formula $C_{29}H_{36}N_2O_6$</p> <p>Mol. weight 508,61 g/mol</p>	 
RL-2120	DBCO-NH ₂	<p>Dibenzocyclooctyne-amine</p> <p>CAS-No. 1255942-06-3</p> <p>Formula $C_{18}H_{16}N_2O$</p> <p>Mol. weight 276,33 g/mol</p>	 
PEG6830	DBCO-dPEG(12)-(5)6-carboxyfluorescein	<p>Dibenzoazacyclooctyne-dodeca(ethylene glycol)-(5)6-carboxyfluorescein</p> <p>Formula $C_{46}H_{39}N_3O_{10}$</p> <p>Mol. weight 1234,34 g/mol</p>	 
PEG6810	Bromoacetamido-PEG(24)-DBCO	<p>Bromoacetamido-24(ethylene glycol)-amido-dibenzoazacyclooctyne</p> <p>Formula $C_{71}H_{118}BrN_3O_{27}$</p> <p>Mol. weight 1525,6 g/mol</p>	 
PEG6800	Bromoacetamido-PEG(12)-DBCO	<p>Bromoacetamido-dodeca(ethylene glycol)-amido-dibenzoazacyclooctyne</p> <p>Formula $C_{47}H_{70}BrN_3O_{17}$</p> <p>Mol. weight 996,97 g/mol</p>	 
PEG6790	Bromoacetamido-PEG(4)-DBCO	<p>Bromoacetamido-tetra(ethylene glycol)-amido-dibenzoazacyclooctyne</p> <p>Formula $C_{31}H_{38}BrN_3O_7$</p> <p>Mol. weight 644,55 g/mol</p>	 
PEG6780	DBCO-PEG(24)-MAL	<p>Dibenzoazacyclooctyne-24(ethylene glycol)-maleimide</p> <p>Formula $C_{76}H_{122}N_4O_{29}$</p> <p>Mol. weight 1555,79 g/mol</p>	 
PEG6770	DBCO-PEG(12)-MAL	<p>Dibenzoazacyclooctyne-dodeca(ethylene glycol)-maleimide</p> <p>CAS-No. 2011777-01-6</p> <p>Formula $C_{52}H_{74}N_4O_{17}$</p> <p>Mol. weight 1027,16 g/mol</p>	 

		Product details	Product details
PEG6760	DBCO-PEG(24)-TFP	<p>Dibenzoazacyclooctyne-24(ethylene glycol)-propionyl 2,3,5,6-tetrafluorophenol ester</p> <p>CAS-No. 2754372-40-0</p> <p>Formula $C_{77}H_{116}F_4N_2O_{28}$</p> <p>Mol. weight 1595,75 g/mol</p>	 
RL-2870	ICG-DBCO	<p>Indocyanine green dibenzoazacyclooctyne</p> <p>Formula $C_{63}H_{64}N_4O_5S$</p> <p>Mol. weight 989,27 g/mol</p>	 
LS-4270	Biotin-DBCO	<p>(3aS,4S,6aR)-N-[3-(11,12-Didehydrodibenz[b,f]azocin-5(6H)-yl)-3-oxopropyl]hexahydro-2-oxo-1H-thieno[3,4-d]imidazole-4-pentanamide</p> <p>CAS-No. 1418217-95-4</p> <p>Formula $C_{28}H_{30}N_4O_3S$</p> <p>Mol. weight 502,63 g/mol</p>	 
RL-3670	Halo-DBCO	<p>N-[2-[2-[(6-chlorohexyl)oxy]ethoxy]ethyl]-gamma-oxo-dibenz[b,f]azocene-5(6H)-butanamide</p> <p>CAS-No. 1808119-16-5</p> <p>Formula $C_{29}H_{35}ClN_2O_4$</p> <p>Mol. weight 511,06 g/mol</p>	 
RL-2450	DBCO-PEG(5)-COOH	<p>Dibenzoazacyclooctyne-penta(ethylene glycol)-propionic acid</p> <p>Formula $C_{32}H_{40}N_2O_9$</p> <p>Mol. weight 596,67 g/mol</p>	 
RL-2460	DBCO-PEG(5)-NHS	<p>Dibenzoazacyclooctyne-penta(ethylene glycol)-propionic acid succinimidyl ester</p> <p>CAS-No. 1378531-80-6</p> <p>Formula $C_{36}H_{43}N_3O_{11}$</p> <p>Mol. weight 693,74 g/mol</p>	 
RL-2480	DBCO-PEG(3)-BisSulfonThiol-Linker	<p>Dibenzoazacyclooctyne-PEG(3)-BisSulfon-Thiol-Linker</p> <p>Formula $C_{59}H_{69}N_3O_{14}S_3$</p> <p>Mol. weight 1140,39 g/mol</p>	 
RL-2520	Biotin-PEG(4)-DBCO	<p>Dibenzoazacyclooctyne-tetra(ethylene glycol)-biotin</p> <p>CAS-No. 1255942-07-4</p> <p>Formula $C_{39}H_{51}N_5O_8S$</p> <p>Mol. weight 749,92 g/mol</p>	 
RL-4020	DBCO-C6-Alkyne	<p>N-(propargylamido)dipoyl-dibenzoazacyclooctyne</p> <p>Formula $C_{24}H_{22}N_2O_2$</p> <p>Mol. weight 370,45 g/mol</p>	 
PEG8120	Biotin-PEG(4)-SS-DBCO	<p>N-(2-((3-(3-(azadibenzocyclooctyn-1-yl)-3-oxopropylamino)-3-oxopropyl)disulfanyl)ethyl)-1-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-am ide</p> <p>Formula $C_{44}H_{60}N_6O_9S_3$</p> <p>Mol. weight 913,18 g/mol</p>	 
PEG8140	Biotin-PEG(4)-Dde-DBCO	<p>N-(15-(4,4-dimethyl-2,6-dioxocyclohexylidene)-19-oxo-19-(azadibenzocyclooctyn-1-yl)-3,6,9,12-tetraoxa-16-azanonadecyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide</p> <p>CAS-No. 1807512-43-1</p> <p>Formula $C_{47}H_{61}N_5O_9S$</p> <p>Mol. weight 872,08 g/mol</p>	 
ADC1520	DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB	<p>dibenzoazacyclooctyne-cyclobutane-1,1-dicarboxamide-citrullyl-(4-aminobenzyl alcohol)</p> <p>CAS-No. 2576471-51-5</p> <p>Formula $C_{37}H_{40}N_6O_6$</p> <p>Mol. weight 664,75 g/mol</p>	 


ADC1530 DBCO-cyclobutane-1,1-dicarboxamide-Cit-PAB-PNP

dibenzazacyclooctyne-cyclobutane-1,1-dicarboxamide-citrallyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2576471-34-4

Formula $C_{44}H_{43}N_7O_{10}$

Mol. weight 829,85 g/mol



Product details


**ADC1620 DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB**

dibenzazacyclooctyne-cyclobutane-1,1-dicarboxamide-alanyl-(4-aminobenzyl alcohol)

CAS-No. 2576471-46-8

Formula $C_{34}H_{34}N_4O_5$

Mol. weight 578,66 g/mol


**ADC1630 DBCO-cyclobutane-1,1-dicarboxamide-Ala-PAB-PNP**

dibenzazacyclooctyne-cyclobutane-1,1-dicarboxamide-alanyl-(4-aminobenzyl)-(4-nitrophenyl)-carbonate

CAS-No. 2576471-43-5

Formula $C_{44}H_{37}N_5O_9$


Mol. weight 743,76 g/mol

**RL-4110 DBCO-Suc-SS-COOH**

CAS-No. 2749426-25-1

Formula $C_{24}H_{24}N_2O_4S_2$


Mol. weight 468,59 g/mol

**RL-4310 DBCO-Suc-CBT**

N1-(2-cyanobenzo[d]thiazol-6-yl)-N4-(3-(11,12-didehydro-5,6-dihydro-dibenzo[b,f]azocin-yl)-3-oxopropyl) succinamide

Formula $C_{30}H_{23}N_5O_3S$


Mol. weight 533,61 g/mol

**PEG7465 Me-PEG(12)-DBCO**


Methyl-12(ethylene glycol)-amido-dibenzazacyclooctyne

Formula $C_{45}H_{68}N_2O_{14}$

Mol. weight 861,04 g/mol

**RL-4010 DBCO-SNAP**Formula $C_{34}H_{31}N_7O_3$

Mol. weight 585,67 g/mol



Product details

**DACN derivatives**


Product details

RL-2730 DACN(Tos2)N,N'-bis(*p*-toluenesulfonyl)-4,8-diazacyclononyne

CAS-No. 1797508-57-6

Formula $C_{21}H_{24}N_2O_4S_2$


Mol. weight 432,56 g/mol

**RL-2710 DACN(Tos,Ns)**N-(*o*-nitrobenzenesulfonyl)-N'-(*p*-toluenesulfonyl)-4,8-diazacyclononyne

CAS-No. 1797508-58-7

Formula $C_{20}H_{21}N_3O_6S_2$


Mol. weight 463,53 g/mol

**RL-2720 DACN(Tos,Suc-OH)**N-succinoyl-N'-(*p*-toluenesulfonyl)-4,8-diazacyclononyne

CAS-No. 2109751-68-8

Formula $C_{18}H_{22}N_2O_5S$


Mol. weight 378,44 g/mol






**RL-2735 DACN(Tos)*HCl**N-(*p*-toluenesulfonyl)-4,8-diazacyclononyne hydrochloride

CAS-No. 2331322-18-8


Formula $C_{14}H_{18}N_2O_2S^*HCl$


Mol. weight 278,37*36,46 g/mol





		Product details
RL-2725	DACN(Tos,Suc-NHS)	<p>N-(succinoyl-NHS ester)-N'-(<i>p</i>-toluenesulfonyl)-4,8-diazacyclononyne</p> <p>CAS-No. 2411082-26-1 Formula C₂₂H₂₅N₃O₅S Mol. weight 475,52 g/mol</p> 
RL-3630	DACN(Tos,Mal)	<p>N-(maleimidobutyryl)-N'-(<i>p</i>-toluenesulfonyl)-4,8-diazacyclononyne</p> <p>CAS-No. 2411082-28-3 Formula C₂₂H₂₅N₃O₅S Mol. weight 443,52 g/mol</p> 
RL-2737	DACN(Tos2,6-OH)	<p>4,8-Bis(<i>p</i>-toluenesulfonyl)-4,8-diazacyclononyne-6-ol</p> <p>CAS-No. 2109751-74-6 Formula C₂₁H₂₄N₂O₅S₂ Mol. weight 448,55 g/mol</p> 
RL-3600	DACN(Ms)*HCl	<p>N-(Mesyl)-4,8-diazacyclononyne hydrochloride</p> <p>CAS-No. 2331322-16-6 Formula C₈H₁₄N₂O₂S*HCl Mol. weight 202,27*36,46 g/mol</p> 
RL-3610	DACN(Ms,Ns)	<p>N-(Mesyl)-N'-(2-nosyl)-4,8-diazacyclononyne</p> <p>CAS-No. 2411082-25-0 Formula C₁₄H₁₇N₃O₆S₂ Mol. weight 387,43 g/mol</p> 


Tetrazine derivatives












		Product details
RL-2390	MeTz-mPEG (10kDa)	<p>alpha-Methyltetrazine-omega-methoxy-poly(ethylene glycol)</p> <p>Mol. weight 10000 Da</p> 

		Product details
RL-2400	MeTz-mPEG (20kDa)	<p>alpha-Methyltetrazine-omega-methoxy-poly(ethylene glycol)</p> <p>Mol. weight 20000 Da</p> 

		Product details
RL-2130	(Me)Tz-benzoic acid	<p>4-(6-methyl-1,2,4,5-tetrazin-3-yl)benzoic acid</p> <p>CAS-No. 1345866-66-1 Formula C₁₀H₈N₄O₂ Mol. weight 216,2 g/mol</p> 

		Product details
RL-2140	(Me)Tz-butanoic acid	<p>4-(6-methyl-1,2,4,5-tetrazin-3-yl)butanoic acid</p> <p>CAS-No. 1923268-81-8 Formula C₇H₁₀N₄O₂ Mol. weight 182,18 g/mol</p> 

		Product details
RL-2240	Bz-Tz-NHS	<p>2,5-dioxopyrrolidin-1-yl 5-(4-(1,2,4,5-tetrazin-3-yl)benzylamino)-5-oxopentanoate</p> <p>CAS-No. 1244040-64-9 Formula C₁₈H₁₈N₆O₅ Mol. weight 398,37 g/mol</p> 

		Product details	Product details
RL-2230	Bz-(Me)Tz-NHS	<p>2,5-dioxopyrrolidin-1-yl 5-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)benzylamino)-5-oxopentanoate</p> <p>CAS-No. 1454558-58-7</p> <p>Mol. weight 412,41 g/mol</p>	 
RL-2590	Tz-benzoyl-TOTA*TFA	<p>Tz-benzoyl-TOTA*TFA</p> <p>Formula C₁₉H₂₈N₆O₄*C₂H₃O₂</p> <p>Mol. weight 404,46*114,02 g/mol</p>	 
RL-2250	Bz-Tz-PEG(5)-NHS	<p>2,5-dioxopyrrolidin-1-yl 1-(4-(1,2,4,5-tetrazin-3-yl)phenyl)-3-oxo-6,9,12,15,18-pentaoxa-2-azaheneicosan-21-oate</p> <p>CAS-No. 1682653-80-0</p> <p>Formula C₂₇H₃₆N₆O₁₀</p> <p>Mol. weight 604,61 g/mol</p>	 
RL-2580	Tz-benzoic acid	<p>4-(1,2,4,5-tetrazin-3-yl)benzoic acid</p> <p>CAS-No. 1345866-65-0</p> <p>Formula C₉H₆N₄O₂</p> <p>Mol. weight 202,17 g/mol</p>	 
RL-2860	ICG-Tz	<p>Indocyanine green tetrazine</p> <p>Formula C₅₄H₅₇N₇O₄S</p> <p>Mol. weight 900,14 g/mol</p>	 
RL-2330	MeTz-PEG(4)-NHS	<p>Methyltetrazine-PEG(4)-propanoyl succinimidyl ester</p> <p>CAS-No. 1802907-92-1</p> <p>Formula C₂₄H₃₁N₅O₉</p> <p>Mol. weight 533,53 g/mol</p>	 
RL-3905	MeTz-PEG(4)-STP	<p>sodium 2,3,5,6-tetrafluoro-4-((1-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)phenoxy)-3,6,9,12-tetraoxapentadecan-15-oyl)oxy)benzenesulfonate</p> <p>Formula C₂₆H₂₇F₄N₄NaO₁₀S</p> <p>Mol. weight 686,56 g/mol</p>	 
LS-4290	Biotin-PEG(4)-MeTz	<p>N-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)benzyl)-1-(5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamido)-3,6,9,12-tetraoxapentadecan-15-amide</p> <p>CAS-No. 1962919-31-8</p> <p>Formula C₃₁H₄₆N₈O₇S</p> <p>Mol. weight 674,82 g/mol</p>	 
LS-4280	Biotin-MeTz	<p>N-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)benzyl)-5-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide</p> <p>CAS-No. 1802883-51-7</p> <p>Formula C₂₀H₂₅N₇O₂S</p> <p>Mol. weight 427,53 g/mol</p>	 
HAA9170	H-L-Lys(MeTz-PhAc)-OH*TFA	<p>N-(2-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)phenyl)acetyl)-L-lysine TFA salt</p> <p>CAS-No. 2578384-82-2 (net)</p> <p>Formula C₁₇H₂₂N₆O₃*CF₃COOH</p> <p>Mol. weight 358,40*114,02 g/mol</p>	 
HAA9470	H-L-Phe(4-MeTz)-OH*TFA	<p>(S)-2-amino-3-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt</p> <p>CAS-No. 1698038-85-5 net</p> <p>Formula C₁₂H₁₃N₅O₂*CF₃COOH</p> <p>Mol. weight 259,27*114,02 g/mol</p>	 
HAA9480	H-L-Phe(4-Azido-PrTz)-OH*TFA	<p>(S)-2-amino-3-(4-(6-(3-azidopropyl)-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt</p> <p>Formula C₁₄H₁₆N₈O₂*CF₃COOH</p> <p>Mol. weight 328,34*114,02 g/mol</p>	 


HAA9490 H-L-Phe(3-MeTz)-OH*TFA

(S)-2-amino-3-(3-(6-methyl-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt

CAS-No. 2036323-75-6 net

Formula C₁₂H₁₃N₅O₂*CF₃COOH

Mol. weight 259,27*114,02 g/mol



Product details


**HAA9500 H-L-Phe(3-BuTz)-OH*TFA**

(S)-2-amino-3-(3-(6-butyl-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt

CAS-No. 2036323-83-6 net

Formula C₁₅H₁₉N₅O₂*CF₃COOH

Mol. weight 301,35*114,02 g/mol


**HAA9510 H-L-Phe(3-iPrTz)-OH*TFA**

(S)-2-amino-3-(3-(6-isopropyl-1,2,4,5-tetrazin-3-yl)phenyl)propanoic acid trifluoroacetic acid salt

CAS-No. 2421119-12-0 net


Formula C₁₄H₁₇N₅O₂*CF₃COOH

Mol. weight 287,32*114,02 g/mol

**RL-2380 MeTz-mPEG (5kDa)**

alpha-Methyltetrazine-omega-methoxy-poly(ethylene glycol)

Mol. weight 5000 Da


**RL-2340 MeTz-PEG(4)-mal**

Methyltetrazine-PEG(4)-maleimide

CAS-No. 1802908-02-6

Formula C₂₄H₃₀N₆O₇

Mol. weight 514,53 g/mol


**RL-2310 MeTz-PEG(4)-COOH**

Methyltetrazine-PEG(4)-acid

CAS-No. 1802907-91-0

Formula C₂₀H₂₈N₄O₇

Mol. weight 436,56 g/mol


**RL-2360 MeTz-Bzl-NH₂*HCl**

Methyltetrazine-benzylamine*HCl

CAS-No. 1596117-29-1

Formula C₁₀H₁₁N₅*HCl

Mol. weight 201,23*36,46 g/mol



Product details


**RL-2320 MeTz-PhAc-NHS**

Methyltetrazine-phenylacetyl succinimidyl ester

CAS-No. 1644644-96-1

Formula C₁₅H₁₃N₅O₄

Mol. weight 327,29 g/mol


**RL-2110 MeTz-Phenyl-TOTA*TFA**

4-(6-methyl-1,2,4,5-tetrazin-3-yl)-N-(4,7,10-trioxadecane-13-amine)benzamide trifluoroacetate salt

CAS-No. 2250433-74-8 (net)

Formula C₂₀H₃₀N₆O₄*CF₃CO₂H

Mol. weight 418,49*114,02 g/mol


**RL-3915 MeTz-PhAc-Sulfo-NHS**

sodium 1-(2-(4-(6-methyl-1,2,4,5-tetrazin-3-yl)phenyl)acetoxy)-2,5-dioxopyrrolidine-3-sulfonate

CAS-No. 1821017-46-2

Formula C₁₅H₁₂N₅NaO₇

Mol. weight 429,34 g/mol


**RL-2370 MeTz-PEG(4)-NH₂*HCl**

Methyltetrazine-PEG(4)-amine HCl salt

CAS-No. 1802908-05-9 net

Formula C₁₇H₂₅N₅O₄*HCl

Mol. weight 363,41*HCl g/mol


**RL-2300 MeTz-PhAcOH**

Methyltetrazine-phenylacetic acid

CAS-No. 1380500-88-8

Formula C₁₁H₁₀N₄O₂


Mol. weight 230,22 g/mol



Product details

TCO derivatives**TCO1060 TCO-NH₂*HCl**


trans-Cyclooctene-amine hydrochloride
 CAS-No. 1800507-94-1
 Formula C₁₂H₂₂N₂O₂*HCl
 Mol. weight 226,32*36,45 g/mol




Product details

**TCO1000 TCO-NHS**


trans-Cyclooctene succinimidyl carbonate
 CAS-No. 1191901-33-3
 Formula C₁₃H₁₇NO₅
 Mol. weight 267,28 g/mol

**TCO1040 TCO-PEG(4)-COOH**


trans-Cyclooctene-PEG(4)-Acid
 CAS-No. 1802913-21-8
 Formula C₂₀H₃₅NO₈
 Mol. weight 417,49 g/mol

**TCO1010 TCO-PEG(4)-NHS**


trans-Cyclooctene-PEG(4)-carboxy succinimidyl ester
 CAS-No. 1621096-79-4
 Formula C₂₄H₃₈N₂O₁₀
 Mol. weight 514,57 g/mol

**TCO1050 TCO-PEG(3)-mal**


trans-Cyclooctene-PEG(3)-maleimide
 CAS-No. 1809356-72-6
 Formula C₂₆H₄₁N₃O₈
 Mol. weight 523,62 g/mol

**TCO1020 TCO-PEG(12)-NHS**


trans-Cyclooctene-PEG(12)-carboxy succinimidyl ester
 CAS-No. 2185016-39-9
 Formula C₄₀H₇₀N₂O₁₈
 Mol. weight 866,99 g/mol

**TCO1070 TCO-PEG(3)-NH₂*HCl**


trans-Cyclooctene-PEG(3)-amine
 CAS-No. 2028288-77-7
 Formula C₁₉H₃₆N₂O₅*HCl
 Mol. weight 372,51*36,46 g/mol

**Norbornene derivatives****RL-2080 Norbornene-NHS**


(Norbornene-2-yl)-N-hydroxysuccinimidylcarbonate
 CAS-No. 1888335-48-5
 Formula C₁₂H₁₃NO₅
 Mol. weight 251,24 g/mol

**RL-2090 Norbornene-methyl-NHS**

(Norbornene-2-methyl)-N-hydroxysuccinimidylcarbonate
 CAS-No. 1986791-87-0
 Formula C₁₃H₁₅NO₅
 Mol. weight 265,26 g/mol

**HAA9235 H-L-Lys(Norbornene-methoxycarbonyl)-OH*HCl**

N-epsilon-(norbornene-methoxycarbonyl)-L-lysine hydrochloride
 CAS-No. 1378916-76-7 net
 Formula C₁₅H₂₄N₂O₄*HCl
 Mol. weight 296,37*36,46 g/mol



References:

- A Stepwise Huisgen Cycloaddition Process: Copper(I)-Catalyzed Regioselective “Ligation” of Azides and Terminal Alkynes; V. V. Rostovtsev, L. G. Green, V. V. Fokin, K. B. Sharpless; *Angew. Chem. Int. Ed.* 2002; **41**: 2596-2599. [https://doi.org/10.1002/1521-3773\(20020715\)41:14<2596::Aid-anie2596>3.0.co;2-4](https://doi.org/10.1002/1521-3773(20020715)41:14<2596::Aid-anie2596>3.0.co;2-4)
- Cu^I-Catalyzed Alkyne-Azide “Click” Cycloadditions from a Mechanistic and Synthetic Perspective; V. D. Bock, H. Hiemstra, J. H. van Maarseveen; *Eur. J. Org. Chem.* 2006; **2006**: 51-68. <https://doi.org/10.1002/ejoc.200500483>
- A strain-promoted [3 + 2] azide-alkyne cycloaddition for covalent modification of biomolecules in living systems; N. J. Agard, J. A. Prescher, C. R. Bertozzi; *J Am Chem Soc* 2004; **126**: 15046-7. <https://doi.org/10.1021/ja044996f>
- Bioconjugation with strained alkenes and alkynes; M. F. Debets, S. S. van Berkel, J. Dommerholt, A. T. Dirks, F. P. Rutjes, F. L. van Delft; *Acc Chem Res* 2011; **44**: 805-15. <https://doi.org/10.1021/ar200059z>
- Strain-Promoted 1,3-Dipolar Cycloaddition of Cycloalkynes and Organic Azides; J. Dommerholt, F. Rutjes, F. L. van Delft; *Top Curr Chem (Cham)* 2016; **374**: 16. <https://doi.org/10.1007/s41061-016-0016-4>
- Heteroatom-embedded Medium-Sized Cycloalkynes: Concise Synthesis, Structural Analysis, and Reactions; R. Ni, N. Mitsuda, T. Kashiwagi, K. Igawa, K. Tomooka; *Angew. Chem. Int. Ed.* 2015; **54**: 1190-1194. <https://doi.org/10.1002/anie.201409910>
- Biomedical applications of tetrazine cycloadditions; N. K. Devaraj, R. Weissleder; *Acc Chem Res* 2011; **44**: 816-27. <https://doi.org/10.1021/ar200037t>
- Inverse electron demand Diels-Alder (iEDDA)-initiated conjugation: a (high) potential click chemistry scheme; A. C. Knall, C. Slagovc; *Chem Soc Rev* 2013; **42**: 5131-42. <https://doi.org/10.1039/c3cs60049a>
- Highly accelerated inverse electron-demand cycloaddition of electron-deficient azides with aliphatic cyclooctynes; J. Dommerholt, O. van Rooijen, A. Borrmann, C. F. Guerra, F. M. Bickelhaupt, F. L. van Delft; *Nat Commun* 2014; **5**: 5378. <https://doi.org/10.1038/ncomms6378>
- Inverse electron demand Diels-Alder (IEDDA) reactions in peptide chemistry; M. Pagel; *J. Pept. Sci.* 2019; **25**: e3141. <https://doi.org/10.1002/psc.3141>
- Tetrazine-Based Cycloadditions: Application to Pretargeted Live Cell Imaging; N. K. Devaraj, R. Weissleder, S. A. Hilderbrand; *Bioconjug. Chem.* 2008; **19**: 2297-2299. <https://doi.org/10.1021/bc8004446>

Notes



Any Questions or Suggestions?



We are there for you – simply choose one of the numerous possibilities to get in touch!

- 📞 +49 (0) 9231 97121-0
- 📠 +49 (0) 9231 97121-99
- ✉️ info@iris-biotech.de
- 🌐 www.iris-biotech.de

Empowering Peptide Innovation